


GoSimple Workcell for BD FACSymphony™ Cell Analyzers



GoSimple Workcell for BD FACSymphony™ Cell Analyzers

BIOSERO »
GO PARTNER



The GoSimple Workcell for the BD FACSymphony™ Cell Analyzers is a compact, integrated lab automation solution for increasing the sample acquisition throughput with walkaway functionality for your BD FACSymphony™ A1 and A5 SE Cell Analyzers. The GoSimple Workcell comes with predefined hardware, a benchtop table setup, and established methods for quick setup and simple operator use.

WHAT IT CAN DO FOR YOU



Deploy Quickly with Pre-Configured Methods and Coordinated Support

Installs on an existing lab bench with minimal time and effort. Pre-programmed methods require no development. Get 12 months of GoCare for support and access to our online training academy.



Extend Walkaway Time to Maximize Instrument Use

Load up to 44 MTP plates to extend the walkaway time for BD FACSymphony — enabling longer, unattended runs and reducing costly instrument idle time.



Increase Operator Efficiency

Offload tedious, repetitive tasks, and leave more time for knowledge work and creative thinking.

NOTE: The GoSimple Workcell is available for the BD FACSymphony™ A1 and A5 SE Cell Analyzers

BIOSERO'S PROVEN PROCESS SUPPORTS INDUSTRY LEADERS

1,000+
Laboratories

Serving 20 of the
20 Top Pharma

Providing GBG-Scheduled
Workcells Since 2012

Supporting Scientists with
Automation Since 2003

Providing Orchestrated
Solutions Since 2018

3,000+ GBG
Installs

What's Included in the GoSimple Workcell for BD FACSymphony™ Cell Analyzer

SUPPORTS ONE BD FACSYMHPHONY A1 OR BD FACSYMHPHONY A5 SE CELL ANALYZER. THE CELL ANALYZER AND AUTOMATION KIT ARE SOLD AND INSTALLED BY BD.

1. Hardware

- 1 1x Robotic Arm
- 2 2x Plate Hotels
- 3 1x Shaker
- 4 1x Regrip Nest
- 5 1x E-stop
- 6 1x Baseplate with System Controller


Optional: Computer

2. Integration Software

- Green Button Go Scheduler Basic
- Preconfigured methods for quick setup and simple operator use
- 21CFR11 Manager Extension

3. Services

- Installation and setup are provided at lab site
- Biosero GoCare Plus Support Subscription for 12 months from date of shipping, which includes access to Biosero Academy



 **GBG Scheduler**

Do More with Automation

GETTING STARTED IS A BREEZE
WITH GUIDED USER INPUT SCREENS



Wizard built with existing Green Button Go tools




No prior or specialized knowledge required to use




Accepts operator input with parametrized fields



Simple entry point into integrations





BUILD YOUR TEAM'S CAPABILITIES WITH BIOSERO'S ONLINE TRAINING ACADEMY

Learn how to modify your pre-configured Green Button Go Methods

Learn about advanced uses of Green Button Go to plan future automation expansions

Learn about how to teach the robotic arm for placing additional instruments


Get access to wiki articles on Frequently Asked Questions (FAQs)



What GoSimple Can Do For You

GOSIMPLE WORKCELL FOR THE BD FACSYMPHONY™ CELL ANALYZERS SUPPORTS A WIDE RANGE OF APPLICATION AREAS

The Biosero GoSimple Workcell for the BD FACSymphony™ Cell Analyzers enables increased sample acquisition throughput across many application areas such as:



**Drug Discovery
& Development**



Immunology




**Immuno-
Oncology**



**Cell & Gene
Therapy**

Workflow of GoSimple Workcell for the BD FACSymphony™ Cell Analyzer

The GoSimple Workcell for the BD FACSymphony™ Cell Analyzer supports the following workflow steps. These steps can be scheduled to run in continuous cycles.



GOSIMPLE WORKCELL FOR BD FACSYMPHONY™ CELL ANALYZERS – COMPATIBILITY



BD FACSYMPHONY™ A5 SE CELL ANALYZER

Flow cytometer with the flexibility of spectral unmixing or compensation workflows with 5 lasers and 48 detectors for maximum coverage of the fluorochrome emission spectrum



BD FACSYMPHONY™ A1 CELL ANALYZER

Benchtop flow cytometer with flexibility to meet a broad spectrum of research needs from small particle research to 16-color immunophenotyping

SUPPORT FOR YOUR AUTOMATION SOLUTION



Biosero GoCare Subscription

Biosero ensures that you have visibility and control at your fingertips to maximize your investment.

Training at Your Own Pace

Self-serve training courses
and certification to learn
at your own pace

Virtual Troubleshooting

For help via email or
phone with challenges
in your automation

24/7 Access to Biosero Portal

View and manage your
users, organizations,
licenses, and subscriptions
any time of day

Easy to Upgrade and Patch

Easy access to
software upgrades,
patches, and
driver updates



Contact Biosero Today to Start Your Automation Journey

Biosero and Green Button Go are trademarks owned and/or registered by Biosero in the U.S. and/or other countries.

Class 1 Laser Product. For Research Use Only. Not for use in diagnostic or therapeutic procedures.

BD, the BD Logo and BD FACSymphony are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved.