



GoSimple™ Workcell for the BD FACSLyric™ Flow Cytometry System



GoSimple™ Workcell for the BD FACSLyric™ Flow Cytometry System

BIOSERO 



The GoSimple Workcell for the BD FACSLyric™ Flow Cytometry System is a compact, integrated lab automation solution for increasing the sample acquisition throughput with walkaway functionality for your BD FACSLyric™ Flow Cytometry System. The GoSimple Workcell comes with predefined hardware, a benchtop table setup, and established methods for quick setup and simple operator use.

WHAT IT CAN DO FOR YOU



Increase Operator Efficiency

Offload tedious, repetitive tasks, and leave more time for knowledge work and creative thinking.



Improve Data Security and Traceability

Ensure high quality results with 21CFR11 Manager compliance module – supporting secure data handling and traceability.



Extend Walkaway Time to Maximize Instrument Use

Load up to 44 MTP plates to extend the walkaway time for BD FACSLyric – enabling longer, unattended runs and reducing costly instrument idle time.



Get Started Quickly, No Prior Automation Experience Required

Installs on an existing lab bench with minimal time and effort. Pre-programmed methods require no development. Get 12 months of GoCare for support and access to our online training academy.

NOTE: The GoSimple Workcell is available for the BD FACSLyric™ Flow Cytometry System.

BIOSERO'S PROVEN PROCESS SUPPORTS INDUSTRY LEADERS

1,000+ Laboratories

Serving 20 of the 20 Top Pharma

Providing GBG-Scheduled Workcells Since 2012

Supporting Scientists with Automation Since 2003

Providing Orchestrated Solutions Since 2018

3,000+ GBG Installs

What's Included in the GoSimple Workcell for the BD FACSLyric™ Flow Cytometry System

SUPPORTS 1X BD FACSLYRIC™ FLOW CYTOMETRY SYSTEM AND AUTOMATION KIT, PURCHASED SEPARATELY FROM BD


1. Hardware

- 1 1x Robotic Arm
- 2 2x Plate Hotels
- 3 1x Shaker
- 4 1x Regrip Nest
- 5 1x E-stop
- 6 1x Baseplate with System Controller

Optional: Computer

2. Integration Software

- Green Button Go Scheduler
- Preconfigured methods for quick setup and simple operator use
- 21CFR11 Manager Extension





3. Services

- Installation and setup are provided at lab site
- Biosero GoCare Plus Support Subscription for 12 months from date of shipping, which includes access to Biosero Academy

Do More with Automation


GETTING STARTED IS A BREEZE
WITH GUIDED USER INPUT SCREENS

-  Wizard built with existing Green Button Go tools
-  "Lock Down" workflows to prevent unauthorized changes
-  No prior or specialized knowledge required to use
-  Accepts operator input with parametrized fields
-  Simple entry point into integrations



BUILD YOUR TEAM'S CAPABILITIES WITH
BIOSERO'S ONLINE TRAINING ACADEMY

- Learn how to modify your pre-configured Green Button Go Methods
- Learn about advanced uses of Green Button Go to plan future automation expansions
- Learn about how to teach the robotic arm for placing additional instruments
- Get access to wiki articles on Frequently Asked Questions (FAQs)



What GoSimple Can Do For You

GOSIMPLE WORKCELL FOR THE BD FACSLYRIC™ FLOW CYTOMETRY SYSTEM SUPPORTS HIGH-THROUGHPUT CLINICAL WORK

BD has partnered with Biosero Inc., a leader in laboratory automation solutions, to offer our biopharmaceutical, CRO, industry, hospital and reference labs the integration of BD instruments with other laboratory platforms for **fully automated workflows**.



Biopharmaceutical



CRO




Industry



Hospital & Reference Labs

Workflow of GoSimple Workcell for the BD FACSLyric™ Flow Cytometry System

The GoSimple Workcell for the BD FACSLyric Flow Cytometry System supports the following workflow steps. These steps can be scheduled to run in continuous cycles.



GOSIMPLE WORKCELL FOR BD FACSLYRIC™ FLOW CYTOMETRY SYSTEM



BD FACSLYRIC™ FLOW CYTOMETRY SYSTEM

The BD FACSLyric™ Flow Cytometry System is a high-performance, highly sensitive flow cytometer that demonstrates exceptional resolution and improved separation to make dim and rare cell populations easier to resolve.

SUPPORT FOR YOUR AUTOMATION SOLUTION



Biosero GoCare Subscription

Biosero ensures that you have visibility and control at your fingertips to maximize your investment.





Contact Biosero Today to Start the Automation Journey

Biosero and Green Button Go are trademarks owned and/or registered by Biosero in the U.S. and/or other countries. BD FACSLyric™ Flow Cytometers are Class 1 Laser Products. The BD FACSLyric™ Flow Cytometer is for In Vitro Diagnostic Use with BD FACSuite™ Clinical Application for up to six colors. The BD FACSLyric™ Flow Cytometer is for Research Use Only with BD FACSuite™ Application for up to 12 colors. Not for use in diagnostic or therapeutic procedures. BD, the BD Logo and BD FACSLyric are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2025 BD. All rights reserved.